










Name and Model Number	Automatic Wrist Blood Pressure Monitor RX813	Automatic Blood Pressure Monitor RX811	Premium Blood Pressure Monitor RX810	Premium Plus Blood Pressure Monitor RX812
Features				
Automatic Inflation	✓	✓	✓	✓
Time and Date	✓	-	✓	✓
Number of Users	1	1	2	2
Memory	90	30	90	120
Averages Last 3 Readings or Averages All Stored Readings	Averages Last 3 Readings	Averages Last 3 Readings	Averages Last 3 Readings	Averages All Stored Readings
Validated by Hypertension Canada	 Recommended by Recommandé par Hypertension Canada Silver   Argent	 Recommended by Recommandé par Hypertension Canada Silver   Argent	 Recommended by Recommandé par Hypertension Canada Silver   Argent	 Recommended by Recommandé par Hypertension Canada Gold   Or
Cuff Size Included	12.5 - 20.5 cm	Wide Range Cuff (M/L) 22 - 42 cm	Wide Range Cuff (M/L) 24 - 43 cm	Wide Range Cuff (M/L) 22 - 42 cm
Batteries Included	✓	✓	✓	✓
Adapter Included	-	✓	✓	✓
Storage Bag Included	✓	✓	✓	✓
Irregular Heartbeat Detection	✓	✓	✓	-
Hypertension Indicator	✓	✓	✓	✓
Blood Pressure Analyzer Software	-	-	-	✓
Patented Pulse Arrhythmia Technology 	-	-	-	✓
Warranty	Lifetime	Lifetime	Lifetime	Lifetime
UPC	77105811935 7	77105811937 1	77105811936 4	77105811912 8



### Blood Pressure Monitors Validated by Hypertension Canada

Hypertension Canada recommends blood pressure monitors that have had their accuracy validated using the British Hypertension Society (BHS) protocol or European Society of Hypertension (ESH).



Helps detect an irregular or erratic heartbeat, medically referred to as an arrhythmia, can be intermittent and therefore difficult for doctors to detect. If the PAD icon appears frequently during measurements consult your doctor.

**Blood Pressure Hotline:**  
**1-866-536-2289**